


# Application News



# How to work with the 2D Strain Analysis package on the EchoPAC workstation



## Content

HOW TO WORK WITH THE 2D STRAIN ANALYSIS PACKAGECONTENT .....	2
GETTING STARTED.....	4
OPTIMIZING IMAGES .....	4
START THE ANALYSIS.....	4
Define your view .....	5
Define your region of interest.....	5
Verify the tracking quality.....	6
THE CURVES.....	7
The different modes .....	8
THE DIFFERENT SCREENS.....	9
The Trace screen.....	9
The single screen .....	9
The quad screen .....	9
THE RESULTS SCREEN .....	11
Adjust your measurements.....	12
EXPORTING YOUR DATA .....	13
EXIT THE ANALYSIS.....	15

### NOTE

This hand out is additional training material.  
For more information please refer to the user manual and/or reference manual.



## Getting started



1. Create an exam
2. Connect the ECG
  - a. Press Physio and change ECG Lead to get the best signal
  - b. Obtain a stable ECG trace

## Optimizing images

1. Sector width
  - o Not too small, the myocardium must be visible during the entire cardiac cycle.
  - o Not too big, this lowers the frame rate
2. Frame rate
  - o Optimal between 40-90 fps
  - o Optimize the frame rate with the rotary knob
3. Use the dual focus if the image quality in the apex is not acceptable
  - o Check the frame rate and increase it accordingly with the rotary.
4. Store loops from all apical views
  - o Apical 4 CH
  - o Apical 2 CH
  - o Apical LAX
  - o It is recommend to acquire all three apical views sequentially in order to get comparable heart rates in all views.

## Start the analysis

Press the Advanced button on the right side menu to get all options.  
Then you will get the button for the Q-Analysis to click on.




On the bottom you will find the button for 2D Strain.



### Define your view

A new window opens where you have to define the view you want to analyze.



### Define your region of interest

On the next screen you will get your loop focused on end-systole. It is moving a little bit back and forth (Yoyo) to better visualize the endocardial border.

The endocardial border has to be defined by the user, following the instructions on the screen where to start (i.e. 4 CH view – begin on basal Septum)




Make a double click to define the end of the definition.

Note:

Apical views are starting from mitralvalve ring and end up at mitralvalve ring on the opposite side


Short Axis views can start everywhere but needs to be a closed circle.  
For more information and examples click on the bottom of the screen on the bar with the Note.

The system will detect the region of interest now automatically.



### Verify the tracking quality

On the next screen the whole loop is analyzed.  
Each wall is separated into 3 segments.



For each wall there is a scoring displayed from the system. This scoring indicates the quality of tracking.

Segment	basSept	midSept	apSept	apLat	midLat	basLat
	V	V	V	V	V	X

Green = good.

Red = not acceptable

You can overrule the system by clicking into the field with the V or X and change the scoring.

Once you can agree on this press approve to come to the next step.




## The curves

Once approved the system will show a screen with three main components.

The analyzed loop with the color-coded parametric information.

An anatomical M-Mode that displays the color over the whole myocardium.

The curves extracted from a sample volume visible in the middle of each segment as a solid dot.



Anatomical M-Mode



The scale for the analyzed modes is automatically set to the normal value for the respective mode.

To adjust the scale to your needs there is a slider in the middle of the menu.



### The different modes

The different modes can be changed in the menu.



Depending on the view there are different possible directions and there are different modes.

Directions of analysis depending on the analysed view.

Longitudinal	all apical views
Transverse	all apical views
Circumferential	Parasternal short axis
Radial	Parasternal short axis

### Analysis modes

V	velocity
SR	strain rate
S	strain
D	displacement
Rot	rotation
Rot R	rotation rate



## The different screens

There are different analysis screens available:



### The Trace screen

That's the first screen that comes up in the analysis.

Upper left     color-coded parametric information in the loop  
Lower left     anatomical M-Mode  
Right side     local trace curves

### The single screen

Only the color-coded loop is displayed here.




### The quad screen

The quad screen is subdivided into the following portions.

Upper left     color-coded loop  
Upper right     the traces  
Upper right     anatomical M-Mode

Upper left     a still frame where you can change between different options  
Click on the title i.e. Peak Systolic Strain. A list will open where you can make your selection.



For each segment the system detects the values and displays them with numbers inside the segment and the according color.

Here two example for the different options.

*Peak detections:*

For each segment the system is detecting the peak value doesn't matter where it is in time of the heart cycle.

*Time detection:*

For each segment the system is detecting the time that is needed to reach the i.e. peak longitudinal strain.


Note:

There are different modes available for apical and parasternal views.



Apical	Parasternal
Cancel	Cancel
Time to Peak Longitudinal Strain	Time to Peak Circumferential Strain
✓ Peak Systolic Strain	Peak Systolic Strain
Wall Motion Rank	Wall Motion Rank
Time to Peak Longitudinal Strain Rate	Time to Peak Circumferential Strain Rate
Peak Systolic Strain Rate	Peak Systolic Strain Rate
Peak Longitudinal Displacement	Post Systolic Index (Circumferential)
Post Systolic Index	Pre-stretch Index
Pre-stretch Index	Time to Peak Radial Strain
Peak Transverse Displacement	Peak Radial Displacement
	✓ Peak Radial Strain
	Post Systolic Index (Radial)

## The results screen



Upper left image      still frame with the ROI borders  
 Lower left image    anatomical M-Mode  
 Right side            you will see the traces for each segment that was approved with good tracking quality in the beginning (segments with the X will be excluded here).

The table below gives some parameters with automatically detected measurements.

		ATTENTION! Values are averages over segments!					
VL	Approve	VL	basept	midsept	apexpt	midlat	apexpt
SL	Approve	Peak S	5.51	3.99	1.88	2.49	5.43
SrL	Approve	Peak E	-9.97	-8.10	-4.23	2.90	-9.54
DL	Approve	Peak A	-4.60	-3.47	-2.41	2.34	3.90
DT	Approve						6.98
							-13.46
							3.02

On the left side the different modes can be selected.  
Following this nomenclature you will find the following.

Last letter direction of analysis

- L = longitudinal
- T = transverse
- C = circumferential
- R = radial

The first letters analysis mode

- V = velocity
- D = displacement
- S = strain
- Sr = strain rate
- Rot = rotation

The different measurements:

- Velocity
  - Peak S = max positive velocity in systole
  - Peak E = detecting the peak of E-wave
  - Peak A = detecting the peak of A-wave
- Strain
  - Peak G = maximum negative strain in the complete heart cycle
  - Peak S = max. negative strain only during systole
  - Peak P = max positive strain value
- Displacement Peak S = maximum value during systole

### Adjust your measurements

It is possible to change the automatically detected measurement points in the result screen. Click on the point to change and move it to the desired position.  
With the next click the measurement point is set and the new values are entered in the table.

## Exporting your data

There are different possibilities to store the analyzed data.  
In the trace screen you have four options to store the data




By pressing one of the options the data will be stored first on the hard drive of the computer.

- Snapshot complete image of the analysis screen – still frame (JPG)
- Store AVI the loop will be stored as running movie (AVI)
- Store Trace all curve details are stored in a text file and can be used for further statistical analysis (TXT)
- Store Full all curve details are stored in a type of excel format (CSV)

In the result screen a still frame of each curve analyze can be stored by pressing the approval button behind each of the parameters.




Every time you choose one of these options a new window will appear where you can add a name to your file and then store it.




Together with the image also a .CSV file is stored. This file contains the measurement results.

Once you have finished all your storage you might want to copy it to a removable media. You can either use a CD/DVD (needs to be formatted before) or an USB memory stick. To move the data from the hard drive to your media press the Files button.




A new window opens where are different actions to take.




- Step 1      Select the type of sort you want to have your files displayed in the window:  
              *By patient*(all related files to this patient are included)  
              *Single files*(every file will be displayed separately)  
              *By 2D strain study date*(sort by the date of your analysis)
- Step 2      Select the files you want to copy from the list
- Step 3      Press the arrow keys to move the files into the lower part
- Step 4      Select the place where you want to copy your data to
- Step 5      Press Run to start the copy operation

## Exit the analysis

Press the exit button to leave the 2D Strain analysis tool.



The system will ask you if you want to store the data.



If you press yes, the system will remember your settings.

Take the loop back from the clipboard – Q-Analysis – 2D Strain. The system will run through the whole processing until the approval image. The only thing to do is to check and approve the settings.