GE Healthcare **Vivid**Club

# Application News



# Restoring images From a Backup media





# Content

RESTORING IMAGES FROM A BACKUP MEDIA	2
General	. 4
Where are the images stored	4
Insert media	5
Recall the images	5

NOTE

This hand out is additional training material. For more information please refer to the user manual and/or reference manual.



# Restore images from a backup media

## General

For all images that are still on the hard drive and therefore are directly accessible the systems shows a small thumbnail image.



All images that are not on the hard drive any longer will show the following placeholder.



#### Where are the images stored

Open a patient

Go to the patients examination list.

Find the name of the backup media under the **Disk** section.

When there is the Local HD: this means the images are still on the systems hard drive and can be directly accessed.

Operid	ld Date DiagCo		Image	Dis	
ADM	02/11/2001		1	Backup_No_01	
ADM	06/09/2005		V.	Backup_No_06	
ADM	06/09/2005		1	Backup_No_07	
ADM	19/02/2008		V.	Local HD:	
ADM	19/02/2008		V	Local HD:	
ADM	19/02/2008		1	Local HD:	
ADM	20/02/2008			Local HD:	



# GE Healthcare **Vivid**Club

## Insert media

Insert the media that was used for the backup of the images.

# **Recall the images**

#### Go to the Image browser.

The system shows the whole list of images sorted by examination date.



To recall the images from the backup media **click** on **one** of the **placeholders** in a specific examination.

Date 02/11/2001			
Date 06/09/2005			

After selecting an image press the **Analyse** button.

A message will come up that this image is not available locally. One of the following options needs to be selected.

- Restore only the Selected (Only the selected image(s) will be restored)
- Restore All Images of the selected exam (All images of this examination will be restored)
- Restore **Current patient**. (Restore all images from all examinations. May require several back up medias)

#### Press **Ok**.

If the media was not inserted before a message will come up to insert the media now. Insert the media and press Ok.

Depending on the amount of data this can take several minutes.

